

SEB1000


SAFETY & COMFORT

The reliable support for manual washing

- Continuous Leak Test
- Complete Traceability
- 7" Colour Touch Screen Display
- Up to 5 Customized Cycles
- Easy & Guided Procedures
- Automatic Detergent Dosing System



medical.cbc-europe.com





TIME OPTIMIZATION

Using the SEB1000 gives users the opportunity to focus on other activities.

QUICK HOSE REPLACEMENT

Quick and easy hose change

MUSCULOSKELETAL DISORDERS

Automatic watering prevents damage caused by repeated and continuous movements

IMMEDIACY

The intuitive color touch display guides the operator through all steps and ensures constant control

BIOHAZARD

Tubing can be sterilized to minimize the risk of cross-contamination

SAFETY FOR ENDOSCOPES

Automatic leak test that keeps the instruments under pressure throughout the cycle to avoid possible failures

TRACEABILITY

Barcode scanner that simplifies the identification of equipment and operators. Reports can be stored on USB and printed for full traceability

SPECIFICATIONS AND TECHNICAL DATA

Model name	SEB1000				
Weight	kg	10	Voltage	V	220 - 240
Dimensions (WxDxH)	mm	230 / 280 / 320	Frequency	Hz	50
Ambient temperature	°C	-10°C - 40°C	Power consumption (max.)	W	40
Relative humidity (max.)	HR%	30 - 80	Power	mA	300



CBC (Europe) GmbH

Fritz-Vomfelde-Str. 6 • 40547 Düsseldorf



medical.cbc-europe.com

info@cbc-europe.com

Technical changes and errors reserved. 082021 v1.0